

# Single-Cell Gene Expression Profiling using TaqMan® Gene Expression Assays with the C<sub>1</sub>™ Single-Cell Auto Prep System

## PURPOSE

This application note evaluates the use of the TaqMan Gene Expression Assays for Specific Target Amplification (STA) with the C<sub>1</sub> Single-Cell Auto Prep System. Please refer to the protocol entitled, *“Using the C<sub>1</sub>™ Single-Cell Auto Prep System to Capture Cells from Cell Culture and Perform Preamplification Using TaqMan® Assays”* (PN 100-6117) for more information.

## INTRODUCTION

When quantifying mRNA, real-time PCR can be performed as either a one-step reaction, where the entire reaction from cDNA synthesis to PCR amplification is performed in a single step, or as a two-step reaction, in which reverse transcription and PCR amplification occur in separate tubes (Wong et al. 2005). This application note compares the C<sub>1</sub> Single-Cell Auto Prep System using a two-step method to a conventional one-step, in-plate reverse transcription (RT) and STA PCR method to evaluate single-cell gene expression using TaqMan assays. The Ambion® Single Cell-to-C<sub>1</sub>™ and the CellsDirect™ protocols were used to perform RT and targeted preamplification to prepare 48 single-cell cDNA libraries for real-time PCR analysis.

Many of the challenges in a one-step single-cell gene expression workflow stem from three important issues:

- 1. Verification:** Inability to individualize and confirm single, viable cells prior to experimentation and analysis
- 2. Reproducibility:** Intra-assay variability stemming from multiple, liquid-handling steps
- 3. Limit of Detection:** Lack of sensitivity with one-step methods to detect low abundant transcripts due to RNA degradation (Wong et al. 2005).

It is widely understood that a two-step method offers significant advantages, including better sensitivity and reproducibility compared to one-step workflow; however, it is significantly more expensive in a tube or plate-based format. The C<sub>1</sub> System is an automated cell preparation solution that offers several technical

and economic advantages. It standardizes the isolation, verification, lysis, and preparation of 96 single, living cells in parallel. It complements cellular enrichment techniques, such as Fluorescence Activated Cell Sorting (FACS) and allows researchers to visually confirm the capture and viability of target cells after enrichment (Figure 1A). The C<sub>1</sub> System eliminates the complicated liquid handling and process of preamplification methods such as STA, minimizing manual manipulation of cell samples and therefore the variability between experiments. Due to the nanoliter reaction scale, the two-step RT-PCR method implemented on the C<sub>1</sub> System offers significant cost saving on a per reaction basis.

The C<sub>1</sub> System utilizes an Integrated Fluidic Circuit (IFC) to enable capture and visualization of single cells using a confocal microscope. In addition, automated scripts enable fluorescent viability/cytotoxicity staining to verify viability upstream of analysis, thereby




Figure 1. C<sub>1</sub> Integrated Fluidic Circuit (IFC) capture site and multi-step reaction architecture. The C<sub>1</sub> IFC contains 96 single-cell capture sites and utilizes multiple compartments to perform multi-step reactions and to mix reagents without manual manipulation. **A.** This panel shows a single BJ fibroblast cell isolated in a capture site on the IFC. **B.** Inset panel displays enlarged image of cell capture site. **C.** Diagram of the IFC architecture downstream of the capture site, illustrating the individual chambers that carry out the step-wise chemistry required for STA. After the cell is lysed the cellular contents flow into the first chamber. Next, RT reagents flow into the first two chambers. Afterward, preamplification reagent flows through the chamber expanding the reaction volume and the reaction is thermal cycled across all five chambers. Finally, harvest reagent pushes the entire contents of all wells out to the harvest inlet. After a simple dilution step, real-time PCR is performed on the BioMark™ HD System.

eliminating variability in the analysis downstream of the experiment. After capture and staining, the C<sub>1</sub> IFC employs multi-step reaction architecture (Figure 1B) to enable fully automated multi-step reaction chemistry. The C<sub>1</sub> Single-Cell Auto Prep System has incorporated automated thermal and pneumatic controls within a closed environment to standardize the addition, mixing, and thermal cycling of all reagents without requiring manual intervention.

## TECHNICAL BACKGROUND

Studies of single-cell eukaryotic gene expression have revealed that transcription occurs in pulses (Chubb et al. 2006; Livak et al. 2013) resulting in significant stochastic variation in transcript quantities from cell-to-cell (Sanchez-Freire et al. 2012). High-throughput, cost-effective methods are required in order to collect enough data from a sufficient number of cells to characterize the variation inherent in single-cell gene expression; performing such a study using conventional real-time PCR in plates would be cost prohibitive due to the large volume of PCR master mix required (Livak et al. 2013). The C<sub>1</sub> Single-Cell Auto Prep System in tandem with the BioMark HD System provides a complete automated workflow to isolate single cells and perform high-throughput real-time PCR, enabling precise assessment of cell-to-cell transcriptional heterogeneity. This workflow enables an automated real-time PCR protocol that provides a cost effective means to evaluate more individual cells and greater sensitivity to detect less abundant transcripts. This two-step protocol was developed specifically for single-cell gene expression analysis.

## EXPERIMENT

### Performance Evaluation of TaqMan Gene Expression Assays on the C<sub>1</sub> System

This experiment evaluates the sensitivity and reproducibility to detect 48 transcripts in the BJ fibroblast cell line over a broad dynamic range using two different RT-STA workflows. The Ambion Single Cell-to-C<sub>t</sub> two-step protocol was carried out on the C<sub>1</sub> PreAmp IFC for cells with an average diameter of 17-25  $\mu\text{m}$  (PN 100-5758) and compared to the CellsDirect one-step protocol for RT and STA performed on cells pre-sorted by FACS in a 96-well plate. The sensitivity of each method was determined by comparing the C<sub>t</sub> values for each transcript using TaqMan Gene Expression assays for transcripts that varied in abundance. TaqMan assays were chosen to interrogate a broad range of transcripts according to previous empirically determined RNA copy number. Analyses comparing the C<sub>t</sub> values for the assays were performed on the single-cell transcriptional profiles to examine the efficacy of each method.

## Materials and Methods

A T-75 flask of BJ fibroblasts were grown to 80% confluence, trypsinized and split into two equal fractions for experimentation. One fraction was used to evaluate the TaqMan Gene Expression Assays with the Single Cell-to-C<sub>t</sub> protocol on the C<sub>1</sub> System. The remaining cells were sorted by FACS into a 96-well PCR plate for analyses with the CellsDirect workflow. For the Single Cell-to-C<sub>t</sub> C<sub>1</sub> protocol, the RT and preamplification steps were performed using 48 pooled TaqMan assays in conjunction with the Ambion Single Cell-to-C<sub>t</sub> kit (Life Technologies, PN 100-4904 B1 and Protocol C, *Fluidigm® Real-Time PCR Analysis Software User Guide*, PN 68000088).

The Single Cell-to-C<sub>t</sub> protocol has been adapted for the C<sub>1</sub> System, whereby the RT and preamplification steps are carried out in succession (on a nanoliter scale) in accordance with the manufacturer's protocol offering greater sensitivity to detect low abundant transcripts. For the C<sub>1</sub> workflow, the C<sub>1</sub> PreAmp IFC was used to isolate individual BJ fibroblast cells on the C<sub>1</sub> System from a prepared suspension of cells (200 cells/ $\mu\text{l}$ ). Single-cell capture and viability were verified on the C<sub>1</sub> IFC using the Live/Dead® Viability/Cytotoxicity Kit for Mammalian Cells (Life Technologies, PN L3224) by on-chip staining and visualization with fluorescence microscopy prior to lysis.

BJ fibroblasts were also sorted by FACS and processed using the CellsDirect protocol with the same 48 pooled TaqMan assays used in the C<sub>1</sub> workflow. The cells were sorted based on viability using the same Live/Dead Viability/Cytotoxicity Kit stated above (Life Technologies, PN L3224) to obtain live cells for analyses. The RT-STA cDNA amplicons from both the C<sub>1</sub> and CellsDirect experiments were diluted according to manufacturer's instructions, and loaded onto a single Gene Expression Dynamic Array™ 96.96 IFC (PN 68000130) for real-time PCR analysis on the BioMark HD System. Additionally, no-template control (NTC) and tube controls (containing 24 BJ fibroblast cells isolated by a dilution series from a count obtained with a hemocytometer or FACS) were performed for both the C<sub>1</sub> and FACS sorted experiments. The comparison of the ensemble averages for each single-cell experiment was compared to their respective tube controls (Figure 2).

## RESULTS

The heat map of  $C_t$  values in Figure 3 illustrates the gradient of transcript levels from low to highly expressed genes from the real-time PCR analysis of the TaqMan assays performed using either the two-step Single Cell-to- $C_t$  or one-step CellsDirect workflows. The real-time PCR result from the cells prepared with the  $C_t$  protocol exhibited a shift in  $C_t$  of three to nine cycles (Figure 4) across a broad range of transcript levels demonstrating greater sensitivity over the CellsDirect workflow performed in plate. In addition, the single cells run on the  $C_t$  System and CellsDirect in tube format captured 48/48 and 34/48 of the transcripts respectively. Of the low abundance transcripts the CellsDirect workflow in-tube format identified 4/16 transcripts, whereas all 16 were identified with the  $C_t$  System. Tube controls were performed whereby the bulk analysis of cells was compared to the ensemble average of the single cell experiments. The correlation between tube (24 cells) and the ensemble single-cell data for Single Cell-to- $C_t$  and CellsDirect methods were  $R^2 = 0.96$  and  $0.85$  respectively (Figure 2).




Figure 2. Ensemble vs. Chip Comparison. The correlation of the ensemble average of single cells vs. tube was determined for each method ( $C_1$  and Single Cell-to- $C_t$  and in-plate and CellsDirect) by comparing the ensemble average of  $C_t$  values of the single cells with a tube control of 24 cells.

## CONCLUSION

TaqMan® Gene Expression Assays demonstrated improved single-cell transcript sensitivity with the  $C_1$  System combined with the two-step Single Cell-to- $C_t$  workflow on an automated platform at the nanoliter scale. The adaptation of TaqMan Assays to a fully integrated microfluidics workflow minimizes reagent usage and enables a cost effective real-time PCR protocol and reproducible analysis of gene expression at the single-cell level. With the combined ability to verify viable cell capture and implement the highly sensitive two-step RT-STA protocol with TaqMan Assays, the  $C_1$  System enables a comprehensive methodology for single-cell gene expression. Moreover, the high-throughput capabilities of the  $C_1$  System combined with the microfluidic-based single-cell gene expression workflow provides a robust platform for generating large data sets for statistically significant analyses with minimal intra-assay variation.





Figure 3. Heat map of  $C_t$  values for  $C_1$ /Single Cell-to- $C_t$  method and in-plate/CellsDirect method over the same 48 TaqMan Assays. The transcripts are sorted from low to highly expressed genes (left to right).



## REFERENCES

Livak, et al. (2013) Methods for qPCR gene expression profiling applied to 1440 lymphoblastoid single cells. *Methods* 59: 71-79.

Sanchez-Freire, et al. (2012) Microfluidic single-cell real-time PCR for comparative analysis of gene expression patterns. *Nature Protocols* 7(5): 829-838

Chubb et al., (2006) Transcriptional pulsing of a developmental gene. *Curr Biol.* 16(10):1018-25.

Wong, et al. (2005) Real-time PCR for mRNA quantitation. *BioTechniques* 39 (1): 75–8.

Pfaffl et al. (2002) Relative expression software tool (REST) for group-wise comparison and statistical analysis of relative expression results in real-time PCR. *Nucleic Acids Res.*;30(9):e36.

---

© 2013 Fluidigm Corporation. All rights reserved. Fluidigm, the Fluidigm logo, BioMark, Dynamic Array and C<sub>1</sub> are trademarks or registered trademarks of Fluidigm Corporation in the U.S and/or other countries. All other trademarks are the property of their respective owners.

For Research Use Only. Not for use in diagnostic procedures.

100-6858 06/2013

### Corporate Headquarters

7000 Shoreline Court, Suite 100  
South San Francisco, CA 94080 USA  
Toll-free: 1.866.FLUIDLINE | Fax: 650.871.7152  
www.fluidigm.com

### Sales

North America | +1 650 266 6170 | info-us@fluidigm.com  
Europe/EMEA | +33 1 60 92 42 40 | info-europe@fluidigm.com  
Japan | +81 3 3662 2150 | info-japan@fluidigm.com  
China (including Hong Kong) | +86 21 3255 8368 | info-china@fluidigm.com  
Asia | +1 650.266.6170 | info-asia@fluidigm.com  
Latin America | +1 650 266 6170 | info-latinamerica@fluidigm.com

