

## **Declaration of Conformity**

# **DSX®** Automated ELISA System

Name and Address of DYNEX Technologies, Inc. Sullyfield Circle

Manufacturer: Chantilly, VA 20151, USA

Name and Address of the DYNEX Technologies, Inc. Yeoman Gate,

Authorized European Yeoman Way, Worthing, West Sussex BN13

Representative: 3QZ, UK

## Conformity

Dynex Technologies Inc. confirms that the DSX has fulfil the applicable obligations imposed by sections 1 to 5 of Annex III and verifies that the device meets the provisions of the In Vitro Diagnostic Medical Devices Directive 98/79/EC

	Name	GMDN	Classification	GHTF
REF		Code		Classification
65100	DSX Automated ELISA		General IVD	Class A
	System ambient		General IVB	Old3371
65200	DSX Automated ELISA		General IVD	Class A
	System with 2 incubators	56676		
65400	DSX Automated ELISA	00070	General IVD	Class A
	System with 4 incubators			
65078-625	DSX-Revelation Software for		Accessory of a	Class A
	DSX™		General IVD	

**DYNEX TECHNOLOGIES Inc.** 14340 Sullyfield Circle Chantilly VA USA Phone: 703.631.7800 Fax: 703.803.1441 Page **1** of **4** 



### Standards Applied

#### Safety & EMC

- Safety Requirements For Electrical Equipment For Measurement, Control, And Laboratory Use -Part 1: General Requirements [UL 61010-1:2012 Ed.3 +R:15Jul2015]
- Safety Requirements For Electrical Equipment For Measurement, Control And Laboratory Use -Part 2-010: Particular Requirements For Laboratory Equipment For The Heating Of Materials [UL 61010-2-010:2015 Ed.3]
- UL 61010-2-101 Issued: 2015/08/14 Ed: 2 Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 2-101: Particular Requirements For In Vitro Diagnostic (IVD) Medical Equipment
- Electromagnetic compatibility: EN 61326-1:2006 with CFR 47, Part 15 Subpart B and ICES-003-4: 2004 for a Class A Device
- IEC 61326-1Issued: 2012/07/10 Ed: 2 Electrical Equipment for Measurement, Control and Laboratory Use EMC Requirements – Part 1: General Requirements
- EN 61326-2-6:2006 Electrical equipment for measurement, control and laboratory use EMC requirements Part 2-6: Particular requirements In vitro diagnostic (IVD) medical equipment
- CSA C22.2#61010-2-081 Issued:2004/07/01 Ed:1 (R2009) Safety requirements for electrical equipment for measurement, control, and laboratory use Part 2-081: Particular requirements for automatic and semiautomatic laboratory equipment for analysis and other purposes
- Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use –
  Part 2-101: Particular Requirements for In Vitro Diagnostic (IVD) Medical Equipment [CSA C22.2#61010-2-101:2015 Ed.2]
- Safety Requirements For Electrical Equipment For Measurement, Control And Laboratory Use Part 2-010: Particular Requirements For Laboratory Equipment For The Heating Of Materials [CSA C22.2#61010-2-010:2015 Ed.3]
   Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use -

Part 2-101: Particular Requirements for In Vitro Diagnostic (IVD) Medical Equipment [CSA C22.2#61010-2-101:2015 Ed.2]

#### **Other Standards**

- Statutory Instruments 2002 No.618 Consumer Protection
- EN 980:2008 Symbols for use in the labelling of medical devices
- EN ISO 13485:2003 Medical devices Quality management systems Requirements for regulatory purposes
- CEN EN ISO 14971:2012 Medical devices Application of risk management to medical devices
- EN ISO 18113-3:2011 In vitro diagnostic medical devices Information supplied by the manufacturer (labelling) -Part 3: In vitro diagnostic instruments for professional use (ISO 18113-3:2009)
- EN 62304:2006 Medical device software Software life-cycle processes
- EN 62366:2008 Medical devices Application of usability engineering to medical devices

Authorized Signatory:

C. Prouso.

Candice Prowse Director of Quality Assurance & Regulatory Affairs Signed at Dynex Technologies Inc. Chantilly, VA On 2018-05-23



DYNEX TECHNOLOGIES Inc. 14340 Sullyfield Circle Chantilly VA USA Phone: 703.631.7800 Fax: 703.803.1441

Page 2 of 4



## **DSX® CERTIFICATE OF COMPLIANCE TO RoHS 2**

Dynex Technologies Inc. certifies that the DSX automated in ELISA analyzer to the best of our knowledge complies with the requirements of Directive 2011/65/EU, on the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment.

The majority of DSX parts do not contain the following chemicals or they are in amounts below the allowable limits as shown in table 1.

Hazardous Substance:	Maximum Concentration:
Lead	1000 ppm
Mercury	1000 ppm
Cadmium	100 ppm
Hexavalent Chromium	1000 ppm
Polybrominated biphenyls	1000 ppm
Polybrominated diphenyl ethers (PBDE)	1000 ppm

The following parts use a RoHS exemptions:

Part Number	Description	Exemption	
23001915	MF55D1215F;RES 12.1M 1%	6C	
24500550	ASSAY FIBER OPTICS	13A 13B	
528300700	JEDEC XYZ1V1.JED U3	6B	
528300800	PMCD160212 FITTING 1/8 BARB PP	6B	
528300900	SML-LX1206GC-TR LED GREEN	6B	
528300901	EXTRUSION 80 X 40 CROSS MEMBER	6B	
528301000	5710-510-10 SPACER BEARING SS	6B	

6B Lead as an alloying element in aluminum containing up to 0.4% lead by weight. 6C Copper alloy containing up to 4% lead by weight. 13A Lead in white glasses used for optical applications 13B Cadmium and lead in filter glasses and glasses used for reflectance standards

**DYNEX TECHNOLOGIES Inc.** 14340 Sullyfield Circle Chantilly VA USA Phone: 703.631.7800 Fax: 703.803.1441



### CHINA RoHS Directive Restrictive Substances Standard SJ/T11364-2014 Table:

	Lead	Mercury	Cadmium	Hexavalent Chromium	Polybrominated Biphenyls	Polybrominated Diphenyl Ethers
	(Pb)	(Hg)	(Cd)	(Cr6)	(PBB)	(PBDE)
Reader module	Х	0	Х	0	0	0
Washer Module	0	0	0	0	0	0
Main Chassis	0	0	0	0	0	0
Casework	0	0	0	0	0	0
Transport Arms	Х	0	0	0	0	0
Incubator Module	0	0	0	0	0	0
Pipette Module	0	0	0	0	0	0

O: indicates that this toxic or hazardous substance contained in all the homogeneous materials for this part, according to EIP-A, EIP-B, EIP-C is below the limit requirement in GB/T 26572

X: indicates that this toxic or hazardous substance contained in all the homogeneous materials for this part, according to EIP-A, EIP-B, EIP-C is above the limit requirement in GB/T 26572



Candice Prowse

Director of Quality Assurance and Regulatory Affairs

2016-10-25

**DYNEX TECHNOLOGIES Inc.** 14340 Sullyfield Circle Chantilly VA USA Phone: 703.631.7800 Fax: 703.803.1441

Page 4 of 4